


Laryngoskope | Laringoscopios  
Laryngoscopes | Laringoscopi

**MC-INTOSH LARYNGOSCOPE**



**BLADE**

Size	A	B	C	Working Length
0	77mm	53mm	14mm	55mm
1	95mm	72mm	43mm	70mm
2	111mm	86mm	44mm	78mm
3	133mm	109mm	44mm	108mm
4	151mm	129mm	44mm	133mm

**Fiber Optic Laryngoscope**

Plastic case, Halogen illumination. Manufactured from AISI 304 Stainless Steel Halogen Illumination 2.5 V. The Bulb is inside the handle in order to simplify the sterilization of the blades. All blades and handles are interchangeable with other leading manufacturer. All sets can be combined as per your request, and delivered in Green plastic case.

Blades Availability:

- \* Matt Finish
- \* Mirror Finish
- \* Sand Blast Finish

Pls, confirm in your order list




**D1-032**  
Laryngoscope Bulb Small



**D1-033**  
Laryngoscope Bulb Large



Laryngoscope Battery Handle for Adult C type 2 Batteries



**Mc-Intosh Laryngoscope Set of Fiber Optic Halogen Illumination**  
Manufactured From: AISI 304 Stainless Steel Material

**Mc-INTOSH FIBER OPTIC LARYNGOSCOPE SET HALOGEN ILLUMINATION**

AVAILABLE SETS WITH BATTERY HANDLE

**D1-029 PEDIATRIC SET**  
Set of 0,1,2  
Blades  
Supplied in Green plastic case

**D1-030 STANDARD SET**  
Set of 1,2,3  
Blades  
Supplied in Green plastic case

**D1-031 ADULT SET**  
Set of 1,2,3,4  
Blades  
Supplied in Green plastic case

