


FIBER OPTIC LARYNGOSCOPE MILLER



BLADE

Size	A	B	C	Working Length
0	75mm	53mm	26mm	55mm
1	111mm	80mm	26mm	80mm
2	151mm	134mm	28mm	130mm
3	188mm	167mm	28mm	170mm
4	203mm	180mm	28mm	180mm



Laryngoscope Battery
Handle for Adult
C type 2 Batteries



Fig. 1



Fig. 2



Fig. 3



Fig. 4

**FIBER OPTIC
LARYNGOSCOPE MILLER**

**AVAILABLE SETS
WITH
BATTERY HANDLE**

D1-034

PEDIATRIC SET

Set of 0,1,2

Blades

Supplied in

Green plastic case

D1-035

STANDARD SET

Set of 1,2,3

Blades

Supplied in

Green plastic case

D1-036

ADULT SET

Set of 1,2,3,4

Blades

Supplied in

Green plastic case